

SHEPHERD™

PERIPHERAL GUIDEWIRE

**AN INTEGRAL
PART OF
PERIPHERAL
INTERVENTION**



CSI®

CARDIOVASCULAR
SYSTEMS, INC.

AN INTEGRAL PART OF PERIPHERAL INTERVENTION.

Guidewires are the most tactile-based tool you use, forging the pathway to treatment. Our peripheral guidewire portfolio is designed to effectively get to and cross the lesion.

4 GRAM SHEPHERD ES

WORKHORSE/EXTRA SUPPORT


COATING - 7 cm Polymer Jacket with Hydrophilic Coating
TAPER LENGTH - 10 cm
RADIOPAQUE - 3 cm Coil



6 GRAM SHEPHERD

WORKHORSE/HIGH TORQUE


COATING - 35 cm Polymer Jacket with Hydrophilic Coating
TAPER LENGTH - 41.5 cm
RADIOPAQUE - 3 cm Coil



12 GRAM SHEPHERD

COMPLEX CROSSING

COATING/COIL LENGTH - 15.5 cm Hydrophilic Coating and 15 cm Coil
TAPER LENGTH - 44 cm
RADIOPAQUE - 3 cm Coil



30 GRAM SHEPHERD

HIGHLY CALCIFIED

COATING/COIL LENGTH - 15.5 cm Hydrophilic Coating and 15 cm Coil
TAPER LENGTH - 33 cm
RADIOPAQUE - 3 cm Coil



SHEPHERD™

PERIPHERAL GUIDEWIRE

ENHANCED HYDROPHILIC COATING

is more lubricious and improves durability, especially during long cases.*

TIP LOAD ESCALATION AND DE-ESCALATION ABILITY

in both workhorse and crossing models.

PROPRIETARY TWISTER CORE WIRE TECHNOLOGY

for enhanced navigation and 1:1 torque transfer and steerability.

IMPROVED PTFE COATINGS

allow catheters to track easily.*

Shepherd™ Peripheral Guidewire Ordering Guide

MODEL NUMBER	PRODUCT NAME	DIAMETER	TIP LOAD SUPPORT	LENGTH
PSP14S-30004	014 Shepherd ES	.014"	4 grams	300 cm
PSP18S-30004	018 Shepherd ES	.018"	4 grams	300 cm
PSP14S-30006	014 Shepherd	.014"	6 grams	300 cm
PSP18S-30006	018 Shepherd	.018"	6 grams	300 cm
PSP14S-30012	014 Shepherd 12	.014"	12 grams	300 cm
PSP18S-30012	018 Shepherd 12	.018"	12 grams	300 cm
PSP14S-30030	014 Shepherd 30	.014"	30 grams	300 cm
PSP18S-30030	018 Shepherd 30	.018"	30 grams	300 cm

*When compared to predecessor Zilient® Peripheral Guidewire

The Shepherd™ Peripheral guidewires are intended to facilitate the placement and exchange of balloon catheters or other interventional devices within the peripheral vasculature during Percutaneous Transluminal Angioplasty (PTA) or other intravascular interventional procedures.

CAUTION: Federal law (USA) restricts this device to sale by, or on the order of, a physician.

For more information, please contact your local CSI representative or call customer service at **1-877-274-0901.**



CARDIOVASCULAR
SYSTEMS, INC.

1225 Old Hwy 8 NW, St. Paul, MN 55112

T: 651.259.1600
877.274.0901

F: 612.677.3355
www.csi360.com

CSI® and Shepherd™ are trademarks of Cardiovascular Systems, Inc.
All other trademarks cited herein are owned by their respective owner.
©2022 Cardiovascular Systems, Inc. EN-7099.B 0722